



Outpatient Surgery Magazine

August 2006

Products for More Efficient GI Procedures New scope technologies, accessories and processing products.



Barrx Medical

HAL0360 System

(888) 66-BARRX

www.barrx.com

Price: Not provided

FYI. Barrx Medical says its HAL0360 System is a safer and more effective therapy for Barrett's Esophagus. To use the system, the surgeon opens the esophagus with a sizing balloon, then inflates a catheter in the area of the Barrett's and activates the HAL0360 energy generator to deliver rapid, controlled, consistent bursts of ablative energy. This energy removes a thinner-than-1mm layer of diseased esophagus, letting healthy tissue replace the ablated tissue in three to four weeks, Barrx says.